

GERDX[®]-System

G-SURG[®] GmbH

Weinbergstrasse 28, 83370 Seeon, Germany

www.g-surg.com – info@g-surg.com

Agenda


1	GERDX®-System	3
2	Indication, procedure & competition	6
	2.1 GERD – Reflux disease	7
	2.2 GIST – Tumor resection	16
3	GERDX®- System: Summary	18

1. A micro-hydraulic platform for tissue approximation, plication and fixation

GERDX®-System at a glance

Patented
micro
hydraulic
technology

- Hydraulic controller for various applications
- Flexible shaft
- Fully steerable
- Precision movement
- “Full view” at any time based on independently moveable endoscope
- Haptic feedback



Distal hydraulic tip with:

- Highest precision
- Highest strength
- For all flexible applications




- Easy handling & short intervention
- Highest patient safety and performance

1. A micro-hydraulic platform for tissue approximation, plication and fixation


One controller for several applications

- Due to its simple handling, the platform can be used after a short training by everyone who is specialized in the area of interventional endoscopy



1. A micro-hydraulic platform for tissue approximation, plication and fixation

Components for GERD, GIST and EWL – Status: in use




Agenda

1	GERDX [®] -System	3
2	Indication, procedure & competition	6
2.1	GERD – Reflux disease	7
2.2	GIST – Tumor resection	16
3	GERDX [®] -System: Summary	18

2.1 GERD – Gastro-esophageal reflux disease

Medical condition

- GERD – ‘Gastroesophageal reflux disease’ is one of the most prevalent diseases in developed countries
- 10-20% of adult populations suffer from at least weekly reflux symptoms¹⁾
- The number of patients who are affected by GERD is increasing worldwide



Heartburn /
chest pain

Higher risk of esophageal cancer from
daily GERD²⁾

Bad breath /
breathing problems

Irritation
of the esophageal tissue

Dysphagia (swallowing)


2.1 GERD – Gastro-esophageal reflux disease

The indication


Introduction of GERD

- GERD is characterized mainly by heartburn and/or regurgitation
- GERD most commonly occurs when the **lower esophageal sphincter (LES)** opens spontaneously or does not close properly
- Left untreated, a persistent state of GERD may lead to **Barrett's esophagus**, a **precancerous condition** which can turn into esophageal cancer

Normal anatomy



Dysfunctional valve



Risk factors

- Obesity
- Hiatus hernia
- Smoking
- Aerated drinks e.g. Coca Cola
- Alcohol
- Coffee
- Fast food
- Spicy food
- Pregnancy
- Certain medications:
 - Asthma drugs
 - Calcium channel blockers
 - Pain killers
 - Sedatives
 - Antidepressants
- etc.

2.1 GERD – Gastro-esophageal reflux disease

Overview of treatments

Proton Pump Inhibitors (PPIs)

- Acid suppression medication
- Treatment might consist in successive long-term drug intake
- PPIs may cause gastro-intestinal and nutritional adverse effects



- ✓ Easy access
- ✓ Sufficient in weak cases
- ✗ Risk of side effects: gastric polyps, pneumonia, nausea, diarrhea, fatigue, headache, vitamin B₁₂ deficiency, acid deficit, and dizziness
- ✗ Might cause addiction
- ✗ Very costly in the long run

Laparoscopic surgeries

- **Fundoplication:** invented by Dr. Nissen and Toupet. The upper part of the stomach is wrapped around the lower end of the esophagus and stitched in place



- ✓ Effectiveness: However more than 50% of operated patients became long-term PPI users 10-15 years after surgery²⁾
- ✗ Invasive
- ✗ Longer hospitalization
- ✗ Risk of side effects: scars, dysphagia/and uncontrolled flatulence, vagus nerve injury, gas bloat syndrome

Endoscopic fundoplication

- **Transoral Full-Thickness-Plications** of the fundus (upper stomach wall) to create a mechanical barrier



- ✓ Incisionless
- ✓ GERD procedures is reversible
- ✓ Short procedure
- ✓ Long-term effectiveness
- ✓ Currently 1 night stationary, potentially ambulant
- ✗ Post procedure: Temporary sore throat, abdominal pain, chest pain
- ✗ Currently not recommended for hiatal hernia > 2cm


27% of PPI users are not satisfied¹⁾

High aversion by patients due to dimension of treatment

New technology allows a patient friendly method and creates a new market

2.1 GERD – Gastro-esophageal reflux disease

Published risks observed with PPI use



- 1) Dharmarajan, TS, et. al., 2008 Mar; 9(3):162-7.
- 2) Eom, CS, et. al., 2011 Feb 22;183(3):310-9.
- 3) Targownik, LE, et. al., 2008 Aug 12; 179(4): 319-26.
- 4) Cahan, MA, et. al. Surg Endosc 2006 Sep; 20(9):1364-7.
- 5) Ho, PM, et. al., JAMA 2009 Mar 4;301(9):937-44.
- 6) Jalving, M, et. al., Aliment Pharmacol Ther. 2006 Nov 1;24(9):1341-8.
- 7) García Rodríguez LA, et. al.;Clin Gastroenterol Hepatol. 2007 Dec;5(12):1418-23.
- 8) Cundy T, Dissanayake A.; Clin Endocrinol (Oxf). 2008 Aug;69(2):338-41.
- 9) Lombardo, L. et. al.; Clin Gastroenterol Hepatol. 2010 Jun;8(6):504-8.
- 10) Lazarus B, et al; JAMA Intern Med. 2016 Feb 1;176(2):238-46.
- 11) Gomm W et al; JAMA Neurol. 2016 Feb 15. doi: 10.1001/jamaneurol.2015. 4791.

2.1 PPI vs. Laparoscopy vs. GERDX®

Comparison

	Proton Pump Inhibitors (PPIs)			Laparoscopy	Endoscopy
	OTC	Prescribed generics	Prescribed brands		
Targeted	✗	✓	✓	(✓)	✓
	Not severe enough for any surgery	These patients should consider an endoscopic intervention		Depending on contra indications	



- Since GERD is a chronic disease instead of an acute illness, it causes **significant costs for the healthcare systems and private households** due to the required long-term disease management
- **Direct costs associated** with the disease include costs of
 - ✗ PPIs
 - ✗ physician office and hospital visits
 - ✗ surgical costs
 - ✗ costs resulting from the disease, such as Barrett's esophagus and esophageal adenocarcinoma
 - ✗ costs of side effects caused by PPI long-term intake

- With an endoscopic intervention
 - ✓ a more invasive laparoscopic surgery and a
 - ✓ potentially lifelong consumption of pharmaceuticals can be avoided or
 - ✓ at least the dose of PPIs can be significantly reduced
- The endoscopic intervention is
 - ✓ very patient friendly short procedure
 - ✓ improves quality-of-life at a
 - ✓ very attractive 'price-performance' ratio

2.1 GERD procedure

Endoplication with GERDX-System

- The GERDX® including an endoscope is inserted at deep sedation or full-anesthesia into the stomach
- Two suture systems are placed for:
 - gastro-gastric plication
 - tight closure of cardia
 - rebuild the flap valve
- Permanent restoration of the barrier to reflux
- The GERDX® enables:
 - full thickness plication
 - exact, intuitive positioning
 - gather large volumes of tissue
- by positioning the sutures with pledges exactly, safe and – as it is a full-thickness-procedure – durable
- The procedure is intuitive, easy to use and completed within approx. 20 minutes



3. GERDX®-System: Summary

Comments on NDO Plicator publications

- **Previous studies were performed with only one suture until 2008** (e.g. the clinical paper – Endoscopic full thickness 5 years data by Pleskow et al). An exception to this is the study done by Renteln et al. which was published in 2007. They used 2 sutures. In the following years 2008 and 2009 Renteln et al published studies with at least 2 placed sutures (publ. #6 & #10). In all subsequent studies, at least 2 sutures were set.
- **Patient selection was significantly less accurate in the 2000s than it is today** (due to poorer measuring instruments such as 24-hour multichannel impedance ph monitoring). In addition, the Hill classification which is used today, allows a better patient selection.
- Khajanchee et al. published a study (publ. #9) in 2009 in which a statistical determination of the **success factors** of a Plicator application was made. As Yashodan pointed out, there was **no corresponding focus on these patients in previous Plicator studies**.
- If the NDP Plicator would be used **today**, taking the above mentioned factors into account (multiple sutures, **improved patient selection** through **better measure procedures** and knowledge of success factors), **this alone should improve the already good results of the Plicator procedure**. *"In conclusion, patient selection is critical in determining who will achieve success."* (Khajanchee publ. #9)
- Furthermore the **GERDX® is characterized by additional constructional advantages** compared to the NDO Plicator™ System incl. **simplified handling, high-precision suturing, single-use**, ... which should lead to a further improvement of the NDO Plicator™ results.